Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/IL05/000304

International filing date: 17 March 2005 (17.03.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US

Number: 60/570,469


Filing date: 13 May 2004 (13.05.2004)

Date of receipt at the International Bureau: 06 June 2005 (06.06.2005)

Remark: Priority document submitted or transmitted to the International Bureau in

compliance with Rule 17.1(a) or (b)





THEREONIAND STRANGS OF WORREON

TO ALL TO WHOM THUSE; PRESENTS SHALL COMES

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

April 29, 2005

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE UNDER 35 USC 111.

APPLICATION NUMBER: 60/570,469

FILING DATE: *May 13, 2004*

PA 1313015

By Authority of the

COMMISSIONER OF PATENTS AND TRADEMARKS

T. LAWRENCE Certifying Officer

U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE nder the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

PROVISIONAL APPLICATION FOR PATENT COVER SHEET This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No. INVENTOR(S) Residence Family Name or Sumame Given Name (first and middle [if any]) (City and either State or Foreign Country) RAMAT GAY ISRAEL 211 ELAN separately numbered sheets attached hereto Additional inventors are being named on the TITLE OF THE INVENTION (500 characters max) FOR THE TREATMENT MID-WRETHEAL SUPPORT DEVICE H FLEXIBLE AND PREVENTION OF LIMINARY CORRESPONDENCE ADDRESS Direct all correspondence to: INCONTINENCE IN FEMALES **Customer Number:** OR ELAN Firm or 211 Individual Name Address Address Zio State GAN KA MAT City Fax Telephone Country ISRA ENCLOSED APPLICATION PARTS (check all that apply) CD(s), Number_ Specification Number of Pages __ Other (specify) ___ Drawing(s) Number of Sheets Application Data Sheet. See 37 CFR 1.76 METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT FILING FEE Applicant claims small entity status. See 37 CFR 1.27. Amount (\$) A check or money order is enclosed to cover the filing fees. \$ 80 The Director is herby authorized to charge filing fees or credit any overpayment to Deposit Account Number: Payment by credit card. Form PTO-2038 is attached. The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government. No. Yes, the name of the U.S. Government agency and the Government contract number are: [Page 1 of 2] 30.404 Date__ Respectfully submitted, REGISTRATION NO. SIGNATURE_ (if appropriate) Docket Number: TYPED or PRINTED NAME +972-TELEPHONE_ USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT

USE URLY FUR HLING A PROVISIONAL APPLICATION FOR PATENT

This collection of information is required by 37 CFR 1.51. The information is required to obtain or retain a benefit by the public which is to fite (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form another suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Application, Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

_	_	ъ.
_	3	-
=	₹	-
>	2	-
≥	•	=
•	3	_

Signature

PTO/SB/17 (10-03) Approved for use through 07/31/2006. OMB 0651-0032

U.S. Patent and Trademark Unites; U.S. DEPARTMENT OF COMMERCE Under the Paperwork Reduction Act of 1895, no persons are required to respond to a collection of information unless it displays a valid OMB control number.							
		Complete if Known					
FEE TRANSMITTAL		Application Number		mber			
		Filing Date			MARGHATE: 2004 May 11 2004		
for FY 2004		First Named Inventor		ventor	Dr Elan Ziv		
Effective 10/01/2003. Patent fees are subject to annual revision.		Examiner Name					
Applicant claims small entity status. See 37 CFR 1.27		Art Unit					
TOTAL AMOUNT OF PAYMENT (\$) 80		Attorney Docket No.		et No.			
		7 deconie			ALCULATION (continued)		
METHOD OF PAYMENT (check all that apply)							
		3. ADDITIONAL FEES Large Entity, Small Entity					
Deposit Account:		e Fee Fee Fee Fee Pescription					
Deposit		de (\$)	Code (\$)	· .	ee Pald	
Account Number	105				harge - late filing fee or oath harge - late provisional filing fee or		
Deposit Account	105			COVE	er sheet		
Name The Director is authorized to: (check all that apply)		3 130			English specification illing a request for ex parte reexamination		
Charge fee(s) indicated below Credit any overpayments		2 2,520 4 920°	1812 2,5 1804 9		uesting publication of SIR prior to		
Charge any additional fee(s) or any underpayment of fee(s)		7 820		Exa	miner action		
Charge fee(s) indicated below, except for the filling fee	180	5 1,840*	1805 1,8	840°Req Exc	uesting publication of SIR after miner action		
to the above-identified deposit account.	125	110	2251	55 Ext	ension for reply within first month	11	
FEE CALCULATION	125	2 420	2252	210 Ext	ension for reply within second month		
1. BASIC FILING FEE Large Entity Small Entity	125	i3 950	2253	475 Ext	ension for reply within third month		
Fee Fee Fee Fee Description Fee Paid Code (\$)	125	54 1,480			ension for reply within fourth month		
1001 770 2001 385 Utility fling fee	125	55 2,010	2255 1	,005 Ext	ension for reply within fifth month		
1002 340 2002 170 Design filing fee	140		2401		tice of Appeal		
1003 530 2003 265 Plant filing fee	140		2402		ng a brief in support of an appeal quest for oral hearing		
1004 770 2004 385 Reissue filing fee	140	03 290 51 1,510	2403 1451 1		ition to institute a public use proceeding		
1005 160 2005 80 Provisional filing fee 80	145		2452	-	ition to revive - unavoidable		
SUBTOTAL (1) (\$) 80		53 1,330	2453		tition to revive - unintentional		
2. EXTRA CLAIM FEES FOR UTILITY AND REISSUE	150	-	2501	665 Uti	lity issue fee (or reissue)		
Extra Claims below Fee Paid	15	02 480	2502	240 De	sign issue fee		
Total Claims 20** = X = X = X	15		2503		ant issue fee		
Claims Mutiple Dependent	14		1460		tilions to the Commissioner		
	18		1807		ocessing fee under 37 CFR 1.17(q)		
Fee Fee Fee Fee Fee Pescription	18		1806		brnission of Information Disclosure Strnt conding each patent assignment per		
Code (\$) Code (\$) 1202 18 2202 9 Claims in excess of 20	80	21 40	8021	pro pro	perty (times number of properties)		
1201 88 2201 43 Independent claims in excess of 3	18	09 770	2809	385 File (37	ing a submission after final rejection 7 CFR 1.129(a))		
1203 290 2203 145 Multiple dependent claim, if not paid	18	10 770	2810	385 Fo	r each additional invention to be		
1204 86 2204 43 ** Reissue Independent claims	١	204 770	2801		emined (37 CFR 1.129(b)) equest for Continued Examination (RCE)		
over original patent 1205 18 2205 9 ** Reissue claims in excess of 20		801 770 802 900			equest for expedited examination		
1205 18 2205 9 Reissue claims in excess of 20 and over original patent			l .	of	a design application		
			pecify) v Basic Fi	lina Fee I	Paid CIETOTAL (2) (C)		
**or number previously paid, if greater, For Reissues, see above Reduced by basic ruling res rate SUBTOTAL (3) (3)							
SUBMITTED BY (Complete (# opplicable)) Registration No. Telephone +97237369261							
Name (Print/Type) Dr Elan Ziv MD			etion No. Y/Agenti		Telephone +9723736926	1	

WARNING: Information on this form may become public. Credit card information should not be included on this form. Provide credit card information and authorization on PTO-2038.

This collection of information is required by 37 CFR 1.17 and 1.27. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 12 minutes to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the Individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADOR! SEND TO: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

Date

March 14, 2004

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

2

A FLEXIBLE MID-URETHRAL SUPPORT DEVICE FOR THE TREATMENT AND PREVENTION OF URINARY INCONTINENCE IN FEMALES

The present invention relates generally to the field of treatment of urinary incontinence in female patients. More specifically, the present invention relates to a disposable device for use in the treatment of urinary incontinence in women. The invention describes a vaginal disposable device which is inserted and removed in a no-self-touch technique, by the patient herself, using a disposable applicator.

Inventor: Dr Elan Ziv, MD OBGYN, Urogynecologist

Background of the Invention

5

10

15

20

25

30

Urinary incontinence is a widespread problem among females. It is estimated that up to 50% of women occasionally leak urine involuntarily, and that approximately 25% of woman will seek medical advice at some point in order to deal with the problem. Stress incontinence, the most common type of urinary incontinence, refers to the involuntary loss of urine resulting from abdominal pressure rise, occurring during exercise, coughing, sneezing, laughing, etc. When stress incontinence occurs, it is usually the result of the abnormal descent of the urethra and bladder neck below the level of the pelvic floor. While many different factors may contribute to the development of stress incontinence, it is most prevalent among women ages 35-65 and those who have had multiple vaginal deliveries. Stress incontinence is both aggravating and unpleasant for women, and it can also be embarrassing. Many women wear sanitary pads or diapers in order to deal with incontinence, though this is not a real solution to the problem and it can be very inconvenient and unreliable. Surgical treatment may involve securing the paraurethal tissues to the periosteum of the pubic bone or the rectus facia in order to elevate the bladder neck above the pelvic floor and thereby distribute pressure equally to the bladder, the bladder neck, and the mid-urethra. Recently, a procedure known as "TVT" ("Tension Free Vaginal Tape") was developed, in which a mesh tape is implanted underneath midurethra, creating a hammock on which the urethra may kink during physical effort.

[&]quot;A flexible midurethral support device for treatment & prevention of female urinary incontinence"

Dr Elan Ziv, MD OBGYN, Urogynecologist

Page 1 of 7

However, surgery is only suitable for severe cases, and the majority of women experiencing incontinence do not need surgical solutions.

One modality of non-surgical treatment involves the use of devices that are inserted into the vagina, either by a medical practitioner or by the woman herself. Most devices are designed to apply pressure against the bladder neck so as to inhibit or completely block the flow of urine through the urethra. A variety of such devices are known in the art. For example, refer to U.S. Patent No. 5,618,256 to Reimer, entitled, "Device for Arrangement in the Vagina for Prevention of Involuntary Urination with Females and an Applicator for use in Insertion of the Device;" U.S. Patent No. 5,785,640 to Kresch, entitled "Method for Treating Female Incontinence;" U.S. Patent No. 4,920,986 to Biswas, entitled, "Urinary Incontinence Device;" U.S. Patent 5,417,226 to Juma, entitled, "Female Anti-Incontinence Device;" U.S. Patent No. 5,386,836 to Biswas, entitled, "Urinary Incontinence Device;" and U.S. Patent No. 5,007,894 to Enhorning, entitled, "Female Incontinence Device."

15

30

5

10

The existing non-surgical incontinence devices suffer from numerous drawbacks:

- A number of devices are constructed so as to completely block the urethra and thus
 they need to be removed or collapsed in order to allow the woman to urinate, an
 inconvenience for the woman wearing the device.
- To overcome this drawback, vaginal devices have been developed having specialized shapes that do not completely block the bladder neck. These devices tend to be large, uncomfortable, and intrusive. They also tend to cause irritation or soreness to the vagina.
- Such devices are expensive to manufacture, and therefore, they are designed to be reusable and/or to remain in the vagina for an extended period of time. Such devices
 are normally made from large bodies of resilient material, such as plastic or hard
 rubber, in order to preserve their functioning for the required amount of time.
 - Most devices known in the art also tend to be difficult or painful to insert and/or remove. In order to correctly inhibit urine flow, the device needs to be properly positioned in the vaginal canal. As stated previously, a doctor may be required to properly position the device.

[&]quot;A flexible midurethral support device for treatment & prevention of female urinary incontinence"

Dr Elan Ziv, MD OBGYN, Urogynecologist

Page 2 of 7

- In cases where a doctor has to insert the device, the device is adapted for remaining in the vagina for a prolonged period of time. When positioned in the vagina for an extended period of the time, the device may cause vaginal infections, necrosis, or bleeding.
- The device may block or inhibit the flow of normal body secretions through the vagina, and may cause inflammation of the vagina and a foul-smelling discharge.
 - In cases where the device is designed to be inserted by the woman herself, the device often has to be removed, cleaned, and then re-inserted after a predetermined number of hours.

10

15

All vaginal devices so far described or marketed have at least one of the limiting features described above. No vaginal device for controlling urinary incontinence has so far been successfully marketed and used by the woman herself. There is a need for a device for controlling involuntary urination that is disposable, easy and comfortable for a woman to use, that works effectively and reliably, and that is completely sanitary and hygienic.

[&]quot;A flexible midurethral support device for treatment & prevention of female urinary incontinence"

Dr Elan Ziv, MD OBGYN, Urogynecologist

Page 3 of 7

The Invention

The present invention provides a device for the treatment of urinary incontinence females. The device of the present invention is adapted to be disposable, worn only for a maximum of 16 hours and then discarded and replaced with a new device (if needed).

- The device of the present invention is simple and easy to use, and is inserted effortlessly in the same user-friendly and familiar manner that a tampon is inserted into the vagina during menstruation. As opposed to large and intrusive devices of the prior art, the device of the present invention is comfortable, and, once inserted, the woman need not think about it again until it is removed.
- When involuntary urination occurs, it is usually the result of the abnormal descent of the bladder neck and the urethra into a low position, away from the intra-abdominal pressure system. This "hypermobility" is the result of some injury to the support mechanism which normally keeps the urethra and the bladder neck in a raised position, along the backside of the pubic bone. The lowering of the bladder neck and the urethra that occur, for example, when a woman coughs, sneezes, or laughs, causing involuntary leakage of urine. The device of the present invention is designed so as to provide a "cradle" or shelf-like support to the urethra whenever the urethra descends momentarily, so as to prevent the leakage of urine. The device does not put pressure against the urethra or the bladder neck, but only provides support when there is a rise in abdominal pressure.

20

25

30

5

The present invention relates to a disposable device for the prevention of involuntary urination in females, adapted for being inserted into the vagina, comprising;

- (a) an internal support structure
- (b) a cover covering said internal support structure and comprised of a flexible material, and;
- (c) an applicator coupled to the internal support structure and the cover for facilitating insertion of the device into the vagina;

The main device and the cover are adapted for forming a cradle support for the midurethra following insertion of the device into the vagina so as to prevent involuntary urination while allowing for voluntary urination.

^{*}A flexible midurethral support device for treatment & prevention of female urinary incontinence"

Dr Elan Ziv, MD OBGYN, Urogynecologist

Page 4 of 7

The invention will now be described with reference to accompanying drawings:

FIG. 1A is a side view of the internal support structure.

FIG. 1B is a top view of the internal support structure

FIG. 1C is a bottom view of the internal support structure

FIG. 2A is a perspective view of the covering

5

10

15

20

25

FIG. 2B is a perspective view of the main device inside the covering.

FIG. 3 is a perspective view of the applicator.

FIG. 4 is a sectional view of the invention within the applicator

FIG.5 is a side view of the female pelvis.

FIG.6 is another embodiment of the device with different shapes and curves of the arms.

The core of the device is a one prolonged embodiment (FIG 1A) which has three distinct parts:

- 1. A top section (12) which serves as the "anchoring" element, for stabilizing the device within the vagina,
- 2. A bottom section (14) which serves as the "supporting" element, generating midurethral support,

Each element of the device (FIG 1A+B+C) has 4 flexible arms. These arms of the anchoring element (12), force the device to remain in situ within the vagina, unable to move inwards or outwards, or to rotate. This occurs as a result of the special tendency of vaginal walls to collapse and form an occluded lumen. The flexible arms of the device cause "tenting" of the walls on top of them with resultant sagging of the walls around the intermediate section, thereby stabilizing the device. The arms of the supporting element (14) cause elevation of the tissues around mid-urethra, acting as a hammock. This hammock supports mid-urethra in a tension free manner, much like the TVT operation. The device is introduced into the cover (FIG 2).

Another embodiment of the device, with different shapes and curves of the arms is shown in FIG 6.

[&]quot;A flexible midurethral support device for treatment & prevention of female urinary incontinence"

Dr Elan Ziv, MD OBGYN, Urogynecologist

Page 5 of 7

The cover (FIG 2A) is made of a flexible smooth mesh material (16) designed as small sack with a string (18). FIG 2B shows the device (20) within the tightly closed mesh cover (22). The cover allows for:

- Reduction of the friction between vagina and the device during insertion & removal.
- 5 Reduction of the friction between the applicator and the device during insertion.

10

15

20

25

30

- Pulling the string causes straightening of the cover, straightening of the vaginal walls, allowing for an easy and smooth removal of the device from the vagina.
- Pulling the string causes the arms to fold slightly towards the midline, thereby reducing its size, allowing for an easy and smooth removal of the device from the vagina.
- The mesh of the cover, being stretched between the arms of the device, also serve as a hammock. In a woman who leaks urine during a stressful event (when abdominal pressure rises during coughing, sneezing, etc.), the urethra sags down but meets the hammock in its mid part. That also causes an elevation of the intra urethral pressure with resultant urinary continence.

The applicator serves for insertion of the device into the vagina (FIG 3), as is done when inserting a regular menstrual tampon. The device is kept within the wider part (26) that is inserted into the vagina. When pushing the plunger (28), the device is pushed through the flower like opening (24), allowing for its immediate action once the applicator is removed from the vagina. The string (32) is visible, protruding out of the opening of the plunger (30).

When the device is still within the applicator (FIG 4), its flexible arms (34) converge towards the midline, allowing for the small dimensions and its insertion via a small diameter applicator. After insertion (FIG 5), the flexible arms of the device gain their preintended tension, enlarge the diameter of the device (46) within the vagina (48), anchoring itself under the bladder (40) between the uterine cervix (36) and the pubic bone (38), supporting mid-urethra (42). The string protrudes out of the vaginal introitus (44), as with the regular menstrual tampon, allowing for removal.

[&]quot;A flexible midurethral support device for treatment & prevention of female urinary incontinence"

Dr Elan Ziv, MD OBGYN, Urogynecologist

Page 6 of 7

The invention has in its basic concept the following features:

- Being a disposable device.
- Insertion of the device is always with an applicator.
- Easy & comfortable insertion and removal.
- Being comfortable to wear.

5

10

20

25

- Being hygiene & odorless
- Being a familiar procedure to most female patients as inserting a menstrual tampon.
- Being inserted by the patient herself, in a no-self-touch technique, with a disposable inserter.
- Being removed by the patient herself, in a no-self-touch technique, with the device collapsing and becoming of small size for painless removal.
- Being of high availability, easy to get everywhere, sold as an Over the Counter (OTC) device.
- Being of low cost.
 - Having complete confidentiality, as with the use of menstrual tampons.
 - Having the ability to be removed instantly when needed.
 - No blockage of vaginal discharge.
 - Wide range of diameters

Alternative embodiments of the invention.


- It may be manufactured in different sizes
- It may be made of many flexible materials, such as silicone, polyurethane, etc.
- It can have more or less than 4 arms.
 - The angle between the arms may be changed.
 - Shape and curves of the arms may vary


"A flexible midurethral support device for treatment & prevention of female urinary incontinence"

Dr Elan Ziv, MD OBGYN, Urogynecologist

Page 7 of 7

"A flexible midurethral support device for treatment & prevention of female urinary incontinence" Dr Elan Ziv, MD, OBGYN, Urogynecologist






"A flexible midurethral support device for treatment & prevention of female urinary incontinence" Dr Elan Ziv, MD, OBGYN, Urogynecologist


"A flexible midurethral support device for treatment & prevention of female urinary incontinence" Dr Elan Ziv, MD, OBGYN, Urogynecologist

"A flexible midurethral support device for treatment & prevention of female urinary incontinence" Dr Elan Ziv, MD, OBGYN, Urogynecologist

4/6



"A flexible midurethral support device for treatment & prevention of female urinary incontinence" Dr Elan Ziv, MD, OBGYN, Urogynecologist



"A flexible midurethral support device for treatment & prevention of female urinary incontinence" Dr Elan Ziv, MD, OBGYN, Urogynecologist